

**Titration Curves Worksheet**

As you work through the steps in the lab procedure, record your experimental values and the results on this worksheet. Use the exact values you record for your data to make later calculations.

**Part B - Titration of KHP with NaOH**

Complete the following table.

**Table A. Titration of KHP**

<b>Concentration of KHP solution</b>	<i>M</i>
<b>Volume of KHP solution titrated</b>	mL
<b>Concentration of NaOH solution</b>	<i>M</i>
<b>Calculated <math>V_{\text{eq}}</math> of NaOH solution</b>	mL

From your titration curve, what is the experimental  $V_{\text{eq}}$  for your KHP titration? Label the  $V_{\text{eq}}$  on each copy of your KHP titration curve. Do not forget to subtract the initial buret reading when determining your  $V_{\text{eq}}$ .

What is the percent difference between your theoretical and experimental equivalence volumes?  
(% Difference =  $((\text{Calculated} - \text{Measured}) \cdot 100) / \text{Calculated}$ .)

Complete the following table.

**Table B. Calculated vs Measured pHs for KHP Titration**

mL of 0.20 M NaOH added	Calculated pH (from prelab)	Measured pH (from titration curve)
0.00		
10.00		
15.00		
20.00		
22.00		

What is the experimental  $pK_a$  value for hydrogen phthalate ( $HP^-$  or  $HC_8H_4O_4^-$ ) that you found at the midpoint of your KHP titration curve? Label the  $pK_a$  on each copy of your KHP titration curve.

The accepted value for the  $pK_a$  of  $HP^-$  is 5.408. How does this compare to your experimental value? What is their percent difference?

Based on your endpoint indicated by the phenolphthalein and the equivalence point determined by the titration curve for your KHP titration using the pH electrode, which of the following are true? (Select all that apply.) (*Note: The order of these options may be different in the WebAssign question.*)

- The phenolphthalein changed color rapidly as the pH changed rapidly.
- The phenolphthalein was not necessary since the pH change indicated the equivalence point.
- The phenolphthalein was necessary; otherwise the equivalence point would have been missed.
- The phenolphthalein changed color rapidly after the pH changed rapidly.
- The endpoint and the equivalence point occurred together.
- The phenolphthalein changed color rapidly before the pH changed rapidly.

### Part C - Titration of $\text{Na}_2\text{CO}_3$ with HCl

Complete the following table.

**Table C. Titration of  $\text{Na}_2\text{CO}_3$**

<b>Concentration of <math>\text{Na}_2\text{CO}_3</math> solution</b>	<i>M</i>
<b>Volume of <math>\text{Na}_2\text{CO}_3</math> solution titrated</b>	mL
<b>Concentration of HCl solution</b>	<i>M</i>
<b>Calculated first <math>V_{\text{eq}}</math> of HCl solution</b>	mL
<b>Calculated second <math>V_{\text{eq}}</math> of HCl solution</b>	mL

From your titration curve, what are the experimental first and second  $V_{\text{eq}}$ 's for your  $\text{Na}_2\text{CO}_3$  titration? Label both  $V_{\text{eq}}$ 's on each copy of your  $\text{Na}_2\text{CO}_3$  titration curve. Do not forget to subtract the initial buret reading when determining your  $V_{\text{eq}}$ 's.

first  $V_{\text{eq}}$

second  $V_{\text{eq}}$

What is the percent difference between your theoretical and experimental equivalence volumes?  
(% Difference = ((Calculated – Measured) · 100)/Calculated.)

% difference first  $V_{\text{eq}}$

% difference second  $V_{\text{eq}}$

Complete the following table.

**Table D. Calculated vs Measured pHs for Na<sub>2</sub>CO<sub>3</sub> Titration**

mL of 0.20 M HCl added	Calculated pH (from prelab)	Measured pH (from titration curve)
0.00		
5.00		
10.00		
15.00		
20.00		
22.00		

What are the experimental  $pK_a$  values for carbonic acid (H<sub>2</sub>CO<sub>3</sub>) and hydrogen carbonate (HCO<sub>3</sub><sup>-</sup>) that you found at the midpoints of your Na<sub>2</sub>CO<sub>3</sub> titration curve? Label the  $pK_a$ 's on each copy of your Na<sub>2</sub>CO<sub>3</sub> titration curve.

$pK_{a1}$

$pK_{a2}$

The accepted values for the  $pK_a$ 's of  $H_2CO_3$  and  $HCO_3^-$  are 6.352 and 10.329, respectively. How do these compare to your experimental values? What are their percent differences?

% difference  $pK_{a1}$

% difference  $pK_{a2}$

Based on your endpoint indicated by the methyl orange and the equivalence points determined by the titration curve for your  $Na_2CO_3$  titration using the pH electrode, which of the following are true? (Select all that apply.) (*Note: The order of these options may be different in the WebAssign question.*)

- The methyl orange was not necessary since the pH change indicated the first equivalence point.
- The endpoint and the first equivalence point occurred together.
- The methyl orange was necessary; otherwise the second equivalence point would have been missed.
- The methyl orange was not necessary since the pH change indicated the second equivalence point.
- The endpoint and the second equivalence point occurred together.
- The methyl orange changed color rapidly as the pH changed rapidly at the first equivalence point.
- The methyl orange changed color rapidly as the pH changed rapidly at the second equivalence point.